Specifications

The resuscitator Ambu Mark IV adult is a reusable hand-operated double-walled self-inflating device. It is used to provide or assist pulmonary resuscitation in adults and children with a body weight >15kg (over 3 years old).

Silicone rubber: inner bag, outer bag, valve disc, supporting strap, mask rim (cuff).
- Polysulphone (PSU): patient valve housing, O-ring, patient
connector, mask dome.
- Polyacetal (POM): inlet valve housing, screw cap (two parts),
connector.
POM or Polypropylene (PP): reservoir housing.Aluminium: nameplate.
Maximum ventilation frequency: depends on the inspiratory volume
used. The expansion rate of the bag must be fact enough to comply with
- The expansion rate of the bag must be fast enough to comply with all recommended frequencies.
- Maximum delivered volume 1300ml.
- Dead volume: <5ml
- Leaks (both ways): non-measurable
- Inspiratory resistance at 50 litres/min: approximately -0.35 kPa (-
$3.5 \text{cmH}_2\text{O}$).
- Inspiratory resistance of reservoir bag at 50 litres/min without O_2
supply: approximately -0.1kPa (-1.0cmH ₂ O)
- Expiratory resistance at 50 litres/min: approximately 0.26 kPa (2.6cmH ₂ O)
- Maximum attainable delivery pressure at 0.2 litre/kPa:
approximately 10kPa (100 cmH ₂ O)
- Application temperature: -18 to 50°C (0 to 122°F)
the elasticity of the bag's outer cover limits the airway pressure to
approximately 7kPa (70cmH2O) when squeezed normally with one
hand, as shown in P.
1 resuscitator, complete, reusable, with patient valve, 1 size-5 mask
1 unit, wrapped in a protective plastic bag; packed in a carton with a
multilingual user and maintenance manual.
Manufacturer's name, product name and reference number, serial
number, and CE mark and reference number of notifying body.
Must be multilingual: English, French and Spanish, others when available.
Mark IV resuscitators should not be exposed to direct sunlight or
elevated temperatures during storage.
Never squeeze the bag too hard prior to storage.
Never store the resuscitator in proximity to products containing oil or grease. Temperature range: -40 to 70°C.

Mark IV resuscitator:
Comprised of 3 parts (A,B,C)
Size: length 275mm, diameter 135mm (including patient valve)
Weight: approximately 415g, resuscitator and patient valve (without mask)

C - Patient valve
B Ballon + Inlet valve A

A) Self inflating bag with inlet valve (A)

- Thin-walled outer cover
- Self-expanding inner bag: volume 1500ml
- Inlet valve housing (ISO 32mm connector)

Extra oxygen can be administered by connecting oxygen to the O_2 inlet valve connector.

- Flange lock and connection, towards patient valve.

- Handle, permits the bag to be held with one hand, for regular inflation.

During the normal one-hand operation of the bag, the elasticity of the outer bag prevents the airway pressure from exceeding 7kPa ($70cmH_2O$), approximately.

B) Patient valve (XANEREAM4AV)

One-way valve made in 3 parts with a single shutter; facilitates the easy removal of foreign matter.

a) Valve housing 21.4 with

- Inspiration connector: 32 mm female for connection 11.3

- patient connector: 22/15 ISO 11.4

b) Valve disc. 21.3The guide pin in the patient valve is a visible indicator of correct functioning of the unit. 21.3c) Expiration connector: 30 mm male (ISO) 11.1/21.2

C) Mask XANEREAM1FM5

Size: 5

Specifications: transparent moulded shell with inflatable rim or cushion.

All components can be disassembled for cleaning, chemicals disinfection and/or steam sterilization at 121 to 134°C (refer to manual for instructions) Manufacturer: Ambu®