

Specifications

Materials:


Medical grade PVC, transparent.

Tube size:


ID 7mm, OD 10mm \pm 0.5mm, bubbled every 90 to 100cm.

The bubbles allow the ends of the tube to be cut on the female or male connector to form either a female or male connector end.

To form a female connector:



To form a male connector:



Unit presentation:

30m reel, clean for medical use, single-use, not sterile.

Packaging:

1 carton, with 30m tubing.

Labelling:

Product name, size, reference number, lot number, expiry date, single-use, manufacturer's name and address, CE marked, and user instructions (English, French and Spanish, others when available), and storage recommendations.

Other requirements:

Must comply with *European Medical Directive 93/42/EE*

- Delivered on a reel, can be cut to the exact length require.

- Vacuum resistant: \pm 559mmHg.

- Single-use.

- Not sterile, can be steam sterilized at 121°C.