

## Specifications

Materials: Dome: polysulphone (PSU), transparent  
Cushioned cuff: rubber (black) or chloroprene (blue)

Description: Mask is made in two parts:



- Dome, lined with a cushioned rim (cuff), with a connector: 15mm ISO.

- Cushioned rim: foam filled (not adjustable or inflatable).

Can be dismantled for cleaning and sterilization (autoclave 121°C).

To be used with RESUSCITATOR AMBU, Mark IV, BABY

Manufacturer: Ambu®

Unit presentation: 1 mask, reusable, in a plastic bag.

Packaging: Protective wrapping with manufacturer's name, product description and reference number; lot number; CE mark and reference number of notifying body; and user and maintenance manual. Must be multilingual: English, French and Spanish, others when available.

Other: Comply with *Medical Device Directive 93/42/EEC*, class IIa

AO mask alternatives:

XANEMASKRN MASK, RENDELL-BAKER, No 0, con. 22mm, neonate, reusable

XANEMASKRP MASK, RENDELL-BAKER, No 1, con. 22mm, baby, reusable