Specifications

Robust, basic multifunctional mechanical operation table with 4 sections for general surgery

Description: Rapidly assembled and dismantled, adjustable on uneven ground,

without special tools.

Mechanically adjustable height.

Mechanical decline and incline, and lateral tilting.

Mechanically adjustable inclinations:

Maximum decline: 35°
Maximum incline: 30°
Maximum lateral tilt: 20°

Crankshaft permits minute adjustment to the table.

Radiotranslucent, permitting X-rays.

Size: Total length: 1.90m

80kg

Overall width: 0.515m
Section widths: 0.475m
Height: 0.75 to 1m

accessories:

Total weight including

Unit presentation: 1 complete table, dismounted, protected in transparent plastic; includes

the required accessories for general, abdominal and perineal surgeries

(refer to accessories list above).

Packaging: Packed in a wooden box, fixed to base for safe transport, includes

assembly instruction: English, French and Spanish.

Labelling: - Designation of item.

- Name and address of supplier (manufacturer).

- Serial number.

Other: Quality complies with *Medical Device Directive* 93/42/ EEC



4 sections:

- Headrest: Length: 0.34m

Incline: 90° to $+60^{\circ}$

- Backrest: Length: 0.55m

Incline: -50° to $+80^{\circ}$

- Seat section, fixed: length: 0.44m

- Leg extension: Length: 0.57m

Incline: -90° to $+30^{\circ}$

Section frames: formed from a plastic sandwich board and anodised aluminium, X-ray permeable, unbreakable, insensitive to climatic conditions, 6mm thickness, and 30x40cm dimension.


Includes 4 antistatic, X-ray permeable cushions, 30mm thickness made from black lectrolite material

- Table base has 4 adjustable jacks on large metallic runners for ease of adjust on uneven ground.

Material: satinized stainless steel, lightweight anodised steel alloy, anticorrosion treatment: phosphatised, coated with anti-chip paint.

NOTE: Solutions with high chlorine concentrations, used for cleaning, can erode the frame's surface and lead to rusting.

Each table is delivered equipped with the following accessories:



| 1 | XSEQTAOPS8R | ROD and infusion stand with 2 hooks | Manufacturer Ref |
|---|-------------|--|------------------|
| 2 | XSEQTAOPS8S | Anaesthesia screen, flexible | TAB005 |
| 3 | XSEQTAOPS8H | Shoulder supports; cylindrical, padded, 1 pair | 45210 |
| 4 | XSEQTAOPS8A | Sliding clamps, to hold accessories, 10 pieces | 45310 |
| 5 | XSEQTAOPS8W | WRISTLETS, padded, 1 pair | 45301 |
| 6 | XSEQTAOPS8T | THIGH REST, padded, 1 pair | 45601 |
| 6 | XSEQTAOPS8L | LEG REST, with straps, 1 pair | 45612 |
| 7 | XSEQTAOPS8Y | INSTRUMENT TRAY, 500x400mm | TAB 281 |
| 8 | XSEQTAOPS8A | ARM REST, padded, for infusion | 45310 |
| | XSEQTAOPS8I | SIDE REST, padded 1 pair | 45401 |
| | XSEQTAOPS8X | CASSETTE HOLDER, for X-ray films | TAB 231 C |